

# Specifications

## OPTICAL OXYGEN FLOW-THROUGH CELLS

### 1 SENSOR SPECIFICATIONS

**Only valid in water/gas (typ. air components) for 2-point calibrated sensors at 20°C, 1013mbar absolute pressure, using default measuring parameters/modes!**

Specifications are valid for oxygen flow-through cells (item no.: **OXFTC**, **OXFTC2**), combined oxygen and temperature flow-through cells (item no.: **TOFTC2**) and flow-through cells with removable oxygen robust probe (item no.: **OXFTCR**).

#### 1.1 Gas Phase: partial pressure pO<sub>2</sub> (hPa), volume percent pV (% O<sub>2</sub> gas)

For a calibrated sensor, the partial oxygen pressure pO<sub>2</sub> in units of hPa (equivalent to mbar) is the fundamental oxygen unit measured by the oxygen meter (in gas and water phase).

Specifications		
<b>Measuring Range</b> Optimum Maximum (not specified)	<b>% O<sub>2</sub> gas</b> 0-50% O <sub>2</sub> 0-100% O <sub>2</sub>	<b>hPa</b> 0-500 hPa 0-1000 hPa
<b>Accuracy *</b> at 1% O <sub>2</sub> /10 hPa at 20% O <sub>2</sub> /200 hPa	±0.02% O <sub>2</sub> ±0.2% O <sub>2</sub>	±0.2 hPa ±2 hPa
<b>Resolution</b> at 1% O <sub>2</sub> /10 hPa at 20% O <sub>2</sub> /200 hPa	0.01% O <sub>2</sub> 0.05% O <sub>2</sub>	0.1 hPa 0.5 hPa
<b>Detection Limit</b>	0.02% O <sub>2</sub>	0.2 hPa

\* The absolute accuracy of full range sensors depends on the calibration mode. For 1-point calibrated sensors these values increase due to a decreasing accuracy. More details on request.

## 1.2 Dissolved Oxygen: % air saturation, $\mu\text{mol/L}$ , $\text{mg/L}$ = $\text{ppm}$ , $\text{mL/L}$

Oxygen dissolved in water can be expressed in % air saturation and in concentration units like  $\mu\text{mol/L}$ ,  $\text{mg/L}$  (ppm), and  $\text{mL/L}$ . For details on calculation of dissolved oxygen units from partial pressure readings (interpolation formula based on temperature, atmospheric pressure and salinity), please see the respective sensor/oxygen meter manuals.

Specifications		
<b>Measuring Range</b> Optimum Maximum (not specified)	<b>% air saturation (a.s.)</b> 0-250% a.s. 0-500% a.s.	<b>mg/L (ppm)</b> 0-22 mg/L 0-44 mg/L
<b>Accuracy *</b> at 5% a.s./0.44 mg/L at 95% a.s./8.8 mg/L	$\pm 0.1\%$ a.s. $\pm 1\%$ a.s.	$\pm 0.01$ mg/L $\pm 0.1$ mg/L
<b>Resolution</b> at 5% a.s./0.44 mg/L at 95% a.s./8.8 mg/L	0.05% a.s. 0.25% a.s.	0.005 mg/L 0.025 mg/L
<b>Detection Limit</b>	0.1% a.s.	0.01 mg/L

\* The absolute accuracy of the full range sensors depends on the calibration mode. For 1-point calibrated sensors these values increase due to a decreasing accuracy. More details on request.

## 1.3 General Characteristics

<b>Response Time (t90) ‡</b> Gas Water	<b>OXFTC/OXFTC2/TOFTC2</b> <1 sec <9 sec	<b>OXFTCR</b> <10 sec <20 sec
<b>Tubing Connectors (Luer-Lock)</b>	ID tubing 1.6 or 2.4 mm (item no. <b>OXFTC</b> ) ID tubing 3.2 or 4.0 mm (item no. <b>OXFTC2, TOFTC2</b> ) ID tubing 3.2 or 4.8 mm (item no. <b>OXFTCR</b> )	
<b>Recommended flow rate for liquids</b>	10-100 mL/min (item no. <b>OXFTC</b> ) 20-500 mL/min (item no. <b>OXFTC2, OXFTCR</b> )	
<b>Temperature Range</b>	specified: 0°C (32°F) to 50°C (122°F) ( <b>OXFTC, OXFTC2, TOFTC2, OXFTCR</b> )	
<b>Minimum Lifetime</b>	10,000,000 data points	
<b>Calibration Modes</b>	1-point and 2-point calibration; obligatory to calibrate in gas ( <i>water</i> ) calibration standards for measurements in gas ( <i>water</i> ) samples	

<b>Application Areas</b>	Laboratory, industry, research. <b>NOT</b> for medical or any safety-critical application. <b>NOT</b> for application in humans. <b>NOT</b> for application in food intended for human consumption.
--------------------------	--

‡ Typical response times for 90% signal change of the oxygen sensor. For liquids: measured for the transition from air into a stirred solution of 1% Na<sub>2</sub>SO<sub>3</sub>

## 2 APPLICABILITY AND CROSS-SENSITIVITY

	<b>Applicability</b>	<b>Cross-Sensitivity</b>	<b>NO Cross-Sensitivity</b>
<b>Water/Aqueous solutions</b>	X		
<b>Gas Phase (typ. air components)</b>	X		
<b>Ethanol<sup>1</sup></b>	short-term only (<1h)		
<b>Methanol<sup>1</sup></b>	short-term only (<1h)		
<b>Isopropanol<sup>1</sup></b>	short-term only (<1h)		
<b>Other organic solvents<sup>2</sup></b>		X	
<b>Chlorine gas (Cl<sub>2</sub>), NO<sub>2</sub> gas, bleach</b>		X	
<b>pH 1-14</b>			X
<b>CO<sub>2</sub></b>			X
<b>CH<sub>4</sub></b>			X
<b>H<sub>2</sub>S</b>			X
<b>Any ionic species</b>			X

<sup>1</sup> Only diluted and after conditioning- contact [info@pyroscience.com](mailto:info@pyroscience.com) for more information;

<sup>2</sup> Includes liquid solvents and solvent vapors

### 3 CLEANING, STERILIZATION, STORAGE

<b>Cleaning</b>	3% H <sub>2</sub> O <sub>2</sub> , Soap solution, short-term Ethanol
<b>Sterilization</b>	short-term 70% Ethanol and 70% Isopropanol
<b>Storage</b>	>3 years in darkness at room temperature

#### Contact

**PyroScience GmbH**  
Hubertusstraße 35  
52064 Aachen  
Deutschland

Tel.: +49 (0)241 5183 2210  
Fax: +49 (0)241 5183 2299  
info@pyroscience.com  
www.pyroscience.com